

LIMS Feature Comparison

Yokogawa LIMS vs Our Pharma-LIMS | 31 May 2026

Category	Yokogawa Feature	Yokogawa	Our LIMS	Our Implementation
Common	User definition screen (workflow)	YES	YES	Role-based users: Admin, Analyst, QA Officer, Lab Manager — configured in Master Data > Users
	Data search	YES	YES	Filter/search across Sample Registration, Work Queue, OOS Investigations, and all master data pages
Test Related	Workflow function	YES	YES	Full lifecycle: Registration !' SRF Sign !' Spec Assignment !' Testing !' Results Review !' QA !' CoA !' Dispatch
	Test requests	YES	YES	Work Queue page lists pending tests per analyst; auto-assigned via Specification Engine on SRF signature
	Issuing test instruction and forms	YES	YES	Form Templates auto-selected on registration; Test Execution renders the active form per method
	Manual input	YES	YES	Enter Results in Test Execution: numeric, text, pass/fail per parameter with OOS/OOT auto-detection
	Automated data import from instru	YES	NO	—
	Batch input	YES	NO	—
	Image capture / file upload	YES	YES	File/image upload in Test Execution with description, stored as evidence per logbook entry
	Computing function (calc formula)	YES	YES	CalcFormulaService evaluates parameter formulas server-side (DataTable.Compute); result shown with &i badge
	Data rounding	YES	YES	DecimalPlaces per parameter; ApplyRounding (MidpointRounding.AwayFromZero) applied server-side on result save
	Automated specification judgment	YES	YES	OOS/OOT detection runs server-side on every result save; auto-transitions to PendingQAReview when OOS closed
Additional Test	Retest function	YES	YES	RetestSampleCommand links retest to original sample; Retest button on Released/Rejected samples with reason m...
	Additional test	YES	NO	—
	Routine test planning	YES	YES	Checkpoint module: Time-Based, OperatorScan, ProcessLog, DispatchEvent trigger modes with configurable intervals
	Stability test function	YES	YES	ICH Q1A programme: protocol registration, pull schedule (T3M–T24M), excursion logging, trend data
Test Support	QC function (Trend / Control Charts)	YES	NO	—
	Reagent & inventory management	YES	YES	Reagents & Standards master data with expiry tracking, lot management, near-expiry alerts on Dashboard
	Analyzer calibration management	YES	YES	Calibration schedule, execution recording (INSERT-only audit), follow-up workflow, e-signature approval
	Connection with analytical instruments	YES	NO	—
	Shipment / Release / CoA issuance	YES	YES	QA release !' CoA PDF generation !' Dispatch QC sign-off !' Delivery Order; e-signature at every step
Reporting	Flexible report formatting	YES	NO	—

Category	Yokogawa Feature	Yokogawa	Our LIMS	Our Implementation
System Management	Master data maintenance (via screen)	YES	YES	Full CRUD: Labs, Instruments, Materials, Test Methods, Parameters, Spec Limits, Form Templates, Users, Reagents
	Master data load (bulk import)	YES	NO	—
	Master data version management	YES	YES	DB-trigger INSERT-only audit on all master data; every change logged with user, timestamp, before/after values
	Supporting tools for batch registration	YES	NO	—
Security & Audit	Password management	YES	YES	BCrypt hashing; §11.300 password re-entry verified independently before every electronic signature
	Access privilege management	YES	YES	JWT role claims; page & action-level access enforced server-side — Analyst cannot approve QA actions
	Screen lock for idle timeout	YES	YES	useIdleLock hook (15 min) + IdleLockOverlay with password re-entry; wired in Layout (21 CFR §11.10(d))
	Saving revision history	YES	YES	INSERT-only audit tables for e-sigs, calibration, checkpoint logs; DB triggers block UPDATE/DELETE on audit rows
External Interface	Receive from other systems	YES	NO	—
	Send to other system	YES	NO	—
Others	Multi language support	YES	NO	—
	Tablet PC support	YES	NO	—

Legend:

YES

 Implemented

NO

 Not implemented

